

DOH FACSCanto II - software parameters and fluorochromes

Laser	Collection Filter (bandpass)	Diva Parameter Name	Fluorochromes
488nm Blue	530/30	530/30 Blue	FITC, Alexa 488, eGFP, CFSE, Dye Cycle Green
	585/42	585/42 Blue	PE, PI, Cy3, Sytox-Orange, Dye Cycle Orange
	670LP	670 Blue	PE-Cy5 (also "tricolour" or "cychrome"), 7AAD, PerCP
	780/60	780/60 Blue	PE-Cy7
633nm Red	660/20	660/20 Red	APC, Alexa 633, Alexa 647, TO-PRO-3, Cy5, DRAQ5, DRAQ7, Dye Cycle Ruby
	780/60	780/60 Red	APC-Cy7, DRAQ7, Alexa Fluor 750
405nm Violet	450/50	450/50 Violet	Alexa 405, CFP, Pacific Blue, DAPI, BV421, Celltrace Violet, Dye Cycle Violet, V450
	510/50	510/50 Violet	Qdot 525, Pacific Orange, Krome Orange, LIVE/DEAD Aqua, V500, BV510

- Only common fluorochromes are listed; please consult FACS Lab staff if your fluorochrome is not present.