

PDI FACS Cantoll configurations

Peter Doherty Institute node				
Laser	BP	Typical fluors	Cantoll b (HTS)	Cantoll d (HTS)
405nm violet	450/40	BV421	◆	◆
	510/50	LD aqua	◆	◆
488nm blue	530/30	FITC	◆	◆
	585/42	PE	◆	◆
	670LP	PerCP, BB700, 7AAD	◆	◆
	780/60	PE-Cy7	◆	◆
633nm red	660/20	APC	◆	◆
	780/60	APC-Cy7	◆	◆

	FILTERS
◆	standard

HTS	High-througput sampler (plate mode)
-----	-------------------------------------